

# YelloPort+plus™ 3mm

Resposable™ Laparoscopic  
3mm Port Access System

YelloPort+plus™ 3mm is a Resposable™ high quality, cost-effective laparoscopic port access system designed for use in ultra-minimally invasive surgery.

The aim of 3mm surgery is to substantially reduce post-operative pain and trauma for the patient, promote faster recovery time and improve 'cosmesis' (minimal scarring) through smaller incisions.



## Key Strengths of YelloPort+plus™ 3mm

- + Resposable™ - reusable Cannulas and Obturators, and disposable sterile, pre-lubricated valves
- + Impressive strength despite reduced size
- + Selection of Cannulas and Obturators available
- + Perfect complement to 3mm Logi™ Range Instruments

For further details about YelloPort+plus™ 3mm please contact:

[sales@surginno.co.uk](mailto:sales@surginno.co.uk)

+44(0)113 230 7597

# YelloPort+plus™ 3mm

Resposable™ Laparoscopic  
3mm Port Access System

## 3mm Obturator Options

Product Name	Product Code	Product Description
YelloPort+plus™ Pyramidal Obturator	YT0305501	55mm
	YT0307501	75mm
	YT0309501	95mm
YelloPort+plus™ Pencil Point Obturator	YT0305503	55mm
	YT0307503	75mm
	YT0309503	95mm

## 3mm Cannulae Options

Product Name	Product Code	Product Description
YelloPort+plus™ Cannula	YC0305500	55mm
	YC0307500	75mm
	YC0309500	95mm
YelloPort+plus™ Cannula + Luer	YC0305501	55mm
	YC0307501	75mm
	YC0309501	95mm
YelloPort+plus™ Cannula + Thread	YC0305510	55mm
	YC0307510	75mm
	YC0309510	95mm
YelloPort+plus™ Cannula+Thread+Luer	YC0305511	55mm
	YC0307511	75mm
	YC0309511	95mm

## 3mm Consumables

Product Name	Product Code	Product Description
YelloPort+plus™ Sterile Valve	YA03VSS01	1 per pack
	YA03VSS02	2 per pack
Luerlock Cap	YALLC10	10 per pack

