




UNIMAX

UNIMAX MEDICAL SYSTEMS INC.



UNIMAX BIPOLAR INSTRUMENTS




www.unimaxmeds.com

UNIMAX BIPOLAR INSTRUMENTS



BIPOLAR CUTTING FORCEPS




Jaw-locking design ensures secure grasp of tissues

Special coated jaws and blade remarkably reduce tissue adhesion

- Multi-functional instrument incorporates grasping, dissecting, transecting, and coagulating functions in one
- Fewer instrument exchanges promote efficiency and reduce procedure time

BIPOLAR DISSECTING FORCEPS



- Large opening of dissecting tip
- Rotation wheel allows the shaft to rotate 360°
- Ergonomic handle increases comfort during operation

BIPOLAR INSTRUMENTS

PRODUCT CODE	DESCRIPTION	PACKAGING
FBS133050	Bipolar cutting forceps 5mm/33cm-with US adaptor	5 pcs/box
FBS233050	Bipolar cutting forceps 5mm/33cm-with EU adaptor- two holes	5 pcs/box
FBS333050	Bipolar cutting forceps 5mm/33cm-with valleylab adaptor- two pins	5 pcs/box
FBS433050	Bipolar cutting forceps 5mm/33cm-with PK adaptor- three pins	5 pcs/box
FBS833050	Bipolar cutting forceps 5mm/33cm- with Erbe adaptor	5 pcs/box
FBD133050	Bipolar dissecting forceps 5mm/33cm-with US adaptor	5 pcs/box
FBD233050	Bipolar dissecting forceps 5mm/33cm-with EU adaptor- two holes	5 pcs/box
FBD333050	Bipolar dissecting forceps 5mm/33cm-with valleylab adaptor- two pins	5 pcs/box
FBD433050	Bipolar dissecting forceps 5mm/33cm-with PK adaptor- three pins	5 pcs/box
FBD833050	Bipolar dissecting forceps 5mm/33cm- with Erbe adaptor	5 pcs/box



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.

UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien City, Taipei, Taiwan
T +886-2-8919-1698 F +886-2-8919-1528 www.unimaxmeds.com