

# EcoSac


























## Semi-Automatic Range

Today's healthcare environment demands greater performance at a lower cost. EcoSac delivers.



*Note: EMP50ECO-R for use with 5mm Trocar*

- Cost effective tissue retrieval system
- Easy deployment
- Mouth or bottom first abdominal entry simplifies sac placement
- Large open mouth to facilitate tissue capture
- Secure drawstring closure

| Product                          | Product Order Code | Trocar Diameter (mm) ID | PRODUCT DIMENSIONS  |             |             | No. Per Box   | Global Trade Item Number (GTIN)   |   |
|----------------------------------|--------------------|-------------------------|---------------------|-------------|-------------|---------------|---|---|
|                                  |                    |                         | Mouth Diameter (mm) | Length (mm) | Volume (ml) |               | Code  | Barcode   |
| EcoSacs<br>Semi Automatic Range  | EMP50ECO-R         | 5                       | 38                  | 90          | 50          | 10            | 5060358120008   |    |
|                                  | EMP70ECO-10        | 10                      | 57                  | 100         | 100         | 10            | 5060358120015   |    |
|                                  | EMP70ECO-S         | 10                      | 57                  | 100         | 100         | 10            | 5060358120022   |    |
|                                  | ECO60GREEN-10      | 10                      | 45                  | 100         | 80          | 10            | 5060358121241   |    |
|                                  | ECO61GREEN-10      | 10                      | 45                  | 100         | 80          | 10            | 5060358121258   |    |
|                                  | ECO70GREEN-10      | 10                      | 70                  | 100         | 100         | 10            | 5060358121104   |    |
|                                  | EMP95ECO-10        | 10                      | 70                  | 155         | 190         | 10            | 5060358120244   |    |
|                                  | EMP100ECO-10       | 10                      | 70                  | 144         | 200         | 10            | 5060358120046   |    |
|                                  | EMP100ECO-S        | 10                      | 70                  | 144         | 200         | 10            | 5060358120053   |    |
|                                  | ECO100GREEN-10     | 10                      | 70                  | 144         | 200         | 10            | 5060358121128   |    |
|                                  | EMP110ECO-R-10     | 5                       | 55                  | 135         | 100         | 10            | 5060358120299   |    |
|                                  | EMP140ECO-10       | 12                      | 102                 | 175         | 550         | 10            | 5060358120077   |    |
|                                  | EMP140ECO-S        | 12                      | 102                 | 175         | 550         | 10            | 5060358120084   |    |
|                                  | EMP145ECO          | 12                      | 102                 | 235         | 750         | 5             | 5060358120114   |    |
|                                  | EMP150ECO-5        | 12                      | 102                 | 270         | 750         | 5             | 5060358120121   |    |
|                                  | EMP160ECO-S        | 12                      | 120                 | 240         | 1000        | 5             | 5060358120138   |    |
|                                  | EMP170ECO          | 12                      | 170                 | 170         | 1200        | 5             | 5060358120145   |    |
|                                  | EMP175ECO-W        | 12                      | 115                 | 255         | 1900        | 5             | 5060358120169   |    |
|                                  | EMP180ECO-5        | 12                      | 170                 | 240         | 2000        | 5             | 5060358120176   |    |
|                                  | EMP200ECO          | 12                      | 170                 | 300         | 2500        | 5             | 5060358120190   |   |
| EMP230ECO                        | 12                 | 153                     | 340                 | 3100        | 5           | 5060358120206 |  |   |
| EMP600ECO-5                      | 15                 | 320                     | 410                 | 6000        | 5           | 5060358120282 |  |   |
| Morcellation Containment Systems | EMP400ECO-TMF-6    | 15                      | 320                 | 410         | 4000        | 6             | 5060358121463   |  |
|                                  | EMP400ECO-TTR-3    | 15                      | 320                 | 410         | 4000        | 3             | 5060358121487   |  |
|                                  | EMP200ECO-TMF-6    | 12                      | 250                 | 325         | 2000        | 6             | 5060358121425   |  |