



Surgical
Innovations

LOGI™

Resposable™ Laparoscopic
Instrumentation



www.surginno.com

Our Logic™ reusable handles are available in two designs, a 'Pistol-Grip' handle and our patented 'Horizontal' handle.


Available with or without ratchet, our 'Pistol-Grip' handle provides a conventional yet ergonomic feel in the surgeon's hand. Logi™Range inserts connect to the handle via

a simple easy-to-use 'ball and socket' connection, providing rapid assembly/disassembly for absolute convenience.

With the option of a centrally mounted, 'three-stage' ratchet mechanism which offers exceptional flexibility and versatility, Logic™ 'Horizontal' handles ergonomic

patented design provides a natural extension to the surgeon's hand.

The electrocautery connection, carefully positioned below the wrist reduces cable interference and snagging for ultimate ease of use.





Logi™Cut
Scissor



Logi™Cut Mini
Scissor



Logi™Dissect
Curved Dissector



Logic™
Pistol-Grip Handle



Logic™
Pistol-Grip Handle with Ratchet



Logi™Clinch
Atraumatic Grasper



Logi™Grasp
Fenestrated Grasper



Logi™Flex
Gastric Band Retriever

Product Information

Product Name	Product Code	Description
Logic™ Horizontal Handle	101 - 41000	Without Ratchet
Logic™ Horizontal Handle	101 - 42000	With Ratchet
Logic™ Pistol-Grip Handle	101 - 43000	Without Ratchet
Logic™ Pistol-Grip Handle	101 - 44000	With Ratchet
Logi™Cut Scissor (16mm blade)	120 - 7001 120 - 7000 120 - 7004	240mm Length (Short) 320mm Length (Standard) 450mm Length (Bariatric)
Logi™Cut Mini Scissor (10mm blade)	120 - 7501 120 - 7500 120 - 7502	240mm Length (Short) 320mm Length (Standard) 450mm Length (Bariatric)
Logi™Dissect Curved Dissector	120 - 7901 120 - 7900 120 - 7902	240mm Length (Short) 320mm Length (Standard) 450mm Length (Bariatric)
Logi™Clinch Atraumatic Grasper	120 - 7701 120 - 7700 120 - 7702	240mm Length (Short) 320mm Length (Standard) 450mm Length (Bariatric)
Logi™Grasp Fenestrated Grasper	120 - 7801 120 - 7800 120 - 7802	240mm Length (Short) 320mm Length (Standard) 450mm Length (Bariatric)
Logi™Flex Gastric Band Retriever	120 - 7200	455mm Length (Standard)

Logi™Range laparoscopic instrument system utilises our innovative resposable™ technology to deliver the cost effectiveness of a reusable instrument with the optimum performance of a disposable. Logi™Range single use inserts are sterile-packed, ready for use and fully disposable, thereby minimising the risk of infection.

Available in a range of configurations, at the fore of endoscopic instrument engineering our hand finished and hand assembled Logi™Range inserts are manufactured 'in-house' to our exacting quality specifications.

Exceeding international standards for sharpness, our scissors and graspers integrate with our reusable handles,

incorporating the latest ergonomic design to enhance tactile feedback for high precision cutting and dissection in every procedure.

Logi™Range.
The high performance, cost-effective solution.



Logic™
Pistol Grip Handle

Logic™
Pistol Grip Handle with Ratchet

Surgical Innovations Ltd
Clayton Wood House
6 Clayton Wood Bank
Leeds
United Kingdom
LS16 6QZ

T. +44 (0) 113 230 7597
F. +44 (0) 113 230 7598

www.surginno.com

For sales enquiries please email
sales@surginno.co.uk

For general enquiries please email
si@surginno.co.uk

Brochure Part Number :
120 – 9005 Issue 1

Product is EC certified and
manufactured under a quality system
certified to ISO 9001 and ISO 13845
and complying with US FDA's Quality
System Regulations.

**Think surgery ...
think Surgical Innovations.**



Part of the
**Surgical Innovations
Group Plc**