

LogiCut™



Resposable® Laparoscopic Scissors


For over 30 years Surgical Innovations have been at the forefront of Minimally Invasive Surgery.

As surgical techniques continue to evolve, the demand for solutions which align to those advancements is key, whilst reducing the **impact on the environment**, as societally we become more aware of the footprint the provision of healthcare generates.

Leading the evolution in its peer group, our **Resposable®** technology is perfectly positioned to meet these challenges and crucially, is accessible to all.

Modular in nature, LogiCut™ combines, without compromise, the key characteristics which clinicians expect from a fully reusable or single use Laparoscopic Scissor, consequently creating a **hybrid solution**.

Designed by a pioneering Laparoscopic Surgeon and manufactured with **handcrafted precision** for over 30 years. The result is an ergonomic construct which delivers precise dissection, exceptional tactile feedback and a seamless user experience for all.





Performance

- Sharp scissors insert ensures **precise dissection** is achieved, equivalent to it's peers
- **Ergonomic** design



Quality

- Resposable® inserts ensure **reduced risk** of compromised insulation over time
- Reduced requirement for use of **complex** decontamination processes and the risk of any **cross contamination** through bioburden



Environmental

- **67%¹ CO2e** reduction versus a fully single use alternative per Laparoscopic Cholecystectomy.
- As part of a wider collection of devices, reduce plastic consumption and the subsequent clinical waste generated by up to **85%¹**



Value

- Ability to reduce costs by up to **17%¹** per Laparoscopic Cholecystectomy.
- Transparent total cost ownership model. No **hidden** servicing, repair, or maintenance costs commonly associated with reusable alternatives

1. Environmental impact and life cycle financial cost of hybrid (reusable/single-use) instruments versus single-use equivalents in laparoscopic cholecystectomy: www.medrxiv.org/content/10.1101/2021.03.10.21253268v2

Ordering Information

5mm: Mini Scissors Profile Insert

Product Code	Product Description
120-7501	5mm LogiCut Mini Scissors Insert 240mm
120-7500	5mm LogiCut Mini Scissors Insert 320mm
120-7502	5mm LogiCut Mini Scissors Insert 450mm

**10mm blade length for all Mini Scissors Inserts*

5mm: Standard Scissor Profile Insert

Product Code	Product Description
120-7001	5mm LogiCut Scissors Insert 240mm
120-7000	5mm LogiCut Scissors Insert 320mm
120-7004	5mm LogiCut Scissors Insert 450mm

**16mm blade length for all Standard Scissors Inserts*



101-43000 Pistol Grip W/O Ratchet



OR





101-41000 Horizontal Handle W/O Ratchet

Leading the Evolution

In Sustainable Surgery.



 Clayton Wood House, 6 Clayton Wood Bank, Leeds, LS16 6QZ, UK

 +44 (0)113 230 7597

 si@surginno.co.uk

 www.surginno.com

 [@SurgicalInnovationsLtd](https://www.linkedin.com/company/SurgicalInnovationsLtd)

 [@SurgInno](https://twitter.com/SurgInno)

 [@SurgicalInnovations](https://www.youtube.com/SurgicalInnovations)