








Master E-Sacs

Fully Automatic Tissue Retrieval System



- ➊ Automatically opening sac
- ➋ Stronger arms reinforce mouth opening
- ➌ Drawstring ensures full closure
- ➍ Integral tail aids removal of sac
- ➎ Accepts tissue sizes up to large spleen, ovary and kidney

Product	Product Order Code	Trocar Diameter (mm) ID	PRODUCT DIMENSIONS			No. Per Box	Global Trade Item Number (GTIN)	
			Mouth Diameter (mm)	Length (mm)	Volume (ml)		Code	Barcode
Master E-Sacs Fully Automated Range	EMP80-WG-6	10	64	105	150	6	5060358120602	
	EMP100-WG-6	10	64	150	200	6	5060358120619	
	EMP120-WG-6	10	85	210	550	6	5060358120633	
	EMP120-WGN-6	10	85	210	550	6	5060358121029	
	EMP150-WG-3	12	110	170	950	3	5060358120671	
	EMP160-WG-3	12	110	260	1500	3	5060358120640	
Robo Sacs	EMP112-WG-6	12	80	150	400	6	5060358120701	