

Logi™Cut 3mm

Resposable™ Laparoscopic 3mm Scissors

Logi™Cut 3mm laparoscopic scissors are a Resposable™ high quality, cost-effective solution for ultra-minimally invasive surgery.

The aim of 3mm surgery is to substantially reduce post-operative pain and trauma for the patient, promote faster recovery time and improve 'cosmesis' (minimal scarring) by making smaller incisions.

Key Strengths of Logi™Cut 3mm

- + Impressive rigidity despite reduced size
- + Resposable™ - reusable handles and single use inserts
- + Double-action scissor blades for easy cutting
- + Insulation up to distal tip of outer tube increasing patient and user safety
- + Better visualisation of tip due to curved nature of scissor
- + Fits our Logic™ reusable handles, including the 'Pistol-Grip' handle and patented 'Horizontal' handle
- + Perfect complement to YelloPort+plus™ 3mm laparoscopic port access system



For further details about Logi™Cut 3mm please contact:

sales@surginno.co.uk

+44(0)113 230 7597


Logi™Cut 3mm

Resposable™ Laparoscopic 3mm Scissors

The entire 3mm Logi™ Range is designed to cut, dissect, manipulate and/or cauterize various tissues during endoscopy/ laparoscopy, general, vascular, gynaecological and thoracic surgical procedures.

| Handle options | | |
|---|--------------|---------------------|
| Product Name | Product Code | Product Description |
| Logi™ Horizontal Handle | 101-41000 | Without ratchet |
| Logi™ Horizontal Handle | 101-42000 | With ratchet |
| Logi™ Pistol-Grip Handle (Standard) | 101-43000 | Without ratchet |
| Logi™ Pistol-Grip Handle (Standard) | 101-44000 | With ratchet |
| Logi™ Small Pistol-Grip Handle | 101-43500 | Without ratchet |

| 3mm scissor options | | |
|---------------------|--------------|--------------|
| Product Name | Product Code | Shaft Length |
| Logi™Cut 3mm | 120-7300 | 320mm |
| Logi™Cut 3mm | 120-7301 | 240mm |



For further details about Logi™Cut 3mm please contact:

sales@surginno.co.uk

+44(0)113 230 7597