

Blue SeromaSet

The optimal single use seroma fluid aspiration system

“For safe and fast removal of seroma fluid”

- Closed system reduces risk of bacterial contamination
- Ergonomic handset for easy control with one hand
- Easier and faster aspiration
- Better emptying due to negative pressure suction
- Patient friendly standard of care
- Fits all needle sizes (needle not included)



Innovative products

Qualified product training

Fast delivery

Clinically trusted quality

World class customer service

Affordable

Blue Surgical
your OR partner

Blue SeromaSet

Seroma formation is the most frequent post-operative complication following breast cancer surgery. SeromaSet has been developed to simplify the emptying of seroma fluid.

SeromaSet is a single use aspiration system consisting of an ergonomic handset which is attached to a negative pressure suction bottle either in 400ml or 600ml. SeromaSet offers a complete closed system, which reduces the risk of bacterial contamination since there is no contact with aspirated fluid by either nurse or patient during the procedure. The aspiration needle is not “left open” during the procedure, which again minimizes the risk of infection inside the breast or at the puncture site.

Due to the effective negative pressure suction bottle, it is possible to make a complete emptying of all seroma fluid in the breast, which can reduce the number of aspirations needed for the patient. SeromaSet is easy to operate by a single staff and with just one hand. The optimal seroma aspiration solution for best patient friendly standard of care. For more information, click www.youtube.com/watch?v=lcLdWMdoBPw.



SeromaSet handset, sold as a set or individually

Available in either 400ml or 600ml bottles

“When caring for post-operative breast cancer patients, it is vital that the process of removing seroma fluid is gentle, safe and quick.”

| REF | Description | Pcs / box |
|---------|---|-----------|
| 1060-10 | SeromaSet, 400ml – handset + 400ml bottle, single use | 10 |
| 1060-20 | SeromaSet, 600ml – handset + 600ml bottle, single use | 20 |
| 1061-10 | SeromaSet, 600ml – handset + 600ml bottle, single use | 10 |
| 1061-20 | SeromaSet, 600ml – handset + 600ml bottle, single use | 20 |
| 1050-10 | SeromaSet handset only, single use | 10 |

Blue Surgical
your OR partner

Exclusive International Distribution by:
Blue Surgical ApS • International Sales Division • Virumvej 64 • DK-2830 Virum • Denmark
Phone: +45 21 19 22 30 • E-mail: info@bluesurgical.com • www.bluesurgical.com
Australian Contact: EMT Healthcare Pty Ltd
email: customerservice@emt.com.au • phone: 02 9979 4111 • www.emt.com.au