



LEVEL 2 Surgical Face Mask (splash resistant)

AUSTRALIAN MADE

(ARTG No: 343201 - Australian Standard AS 4381:2015)



Product Code	Product Description	Quantity
▶ FM-SUR-PK-50	Level 2 Surgical Face Mask, Disposable, 17.5 x 9.5cm	Box/50

Mask Specifications:

Surgical Masks meet AS 4381:2015 Level 2 classification

Bacterial Filtration Efficiency (BFE): ≥ 98%

Differential Pressure (Delta P) mm H₂O/cm²: < 5.0

Resistance to penetration by synthetic blood, minimum pressure in mm Hg for pass result: 120mm Hg