





LEVEL 3 Surgical Face Mask - Type IIR (splash resistant)

TGA Approved - CE Approved - FDA Approved

(ARTG id: 333574 - Bioburden Test EN ISO 11737-1:2018)



Product Code	Product Description	Quantity
FS-001	Level 3 Type IIR Surgical Face Mask, Disposable, 17.5 x 9.5cm	Box/50

Mask Specifications:

Surgical Masks meet EN 14683:2019 Type IIR classification

Bacterial Filtration Efficiency (BFE) ≥ 99%

Resistance to penetration by synthetic blood (fluid resistance) 160mm Hg

Differential Pressure (Delta P) mm H2O/cm2 < 6.0